

購買（英文）規格書

COMMODITY DESCRIPTION

關稅分類番號 H.S. NO.	品目番號 ITEM NO.	品目 및 規格 DESCRIPTION	單位 UNIT	數量 QUANTITY
		내시경보관장	set	1

A. Features:

1. Color touchscreen display for quick identification of endoscope storage time
2. Locking doors protect endoscope inventory
3. Barcode scanner for ease of endoscope and operator identification
4. Convenient transport cassette for endoscope protection and recontamination prevention
5. Validated with incoming instrument-grade compressed air which defines compressed air by particles, humidity and oil: shortens the amount of time needed to fully dry the internal channels of a flexible endoscope compared to traditional hanging methods
6. Ensure your flexible endoscopes and remain safe and dry during storage
7. Drying cabinet continuously forces HEPA filtered air through all endoscope channels, ensuring the endoscope is dry and fluid-free
8. Protect flexible endoscopes from damage: cassette system eliminates free swinging of the distal tip and protects the endoscope during transport and storage.
9. Customizable to fit your facility' s needs: space-conscious design available in various configurations to meet the requirements of any facility.
10. Maximize patient safety: Storage uses a cassette to minimize touching and potential recontamination
11. Data managements ensure endoscope drying cycle compliance: continuous data monitoring provides detailed information on the drying status of all endoscopes.

B. Specifications:

1. Dimensions: 190(H) × 60(W) × 65(D)cm

2. Weight: 170kg
3. Operating Conditions: 10°C to 40°C
4. Storing Conditions: -10°C to 60°C
5. Interfaces: Ethernet RJ45
6. Compressed Air Connection: Pressure (29 - 145psi)
Recommended 5bar
7. Power Supply Voltage: 220V, 60Hz (75W)
8. Display: 6" Color Touch Screen
9. IP Code/Electrical Protection Class: IP20/1

C. Configurations:

1. Storage and Drying System (with Standard Acc'y)	1set
2. Barcode Reader with Cable and Holder	1ea
3. Power Code	1ea
4. Pneumatic Tube with Quick Connector	1ea
5. LAN Cable	1ea
6. Transport Cassette	8ea
7. Hook-up	8ea
8. Service Manual (incl. PDF File)	1ea
9. Operating Manual	1ea

D. Remarks:

1. The equipment should be installation by contractor.
2. Three years warranty after the installation and correct operation.
3. Manufactured products must be delivered within one year based on the inspection date.
4. Products that are not scheduled to be discontinued within one year after delivery must be delivered.
5. A report must be summited after conducting a preventive inspection in accordance with the hospital medical equipment inspection regulations within the warranty period.
6. 보증기간 제외품목: 없음